

Panasonic

RELIABLE CLEANING FOR A RADIANT SMILE.

Our knowledge for your well-being.
The DentaCare oral hygiene system
from Panasonic.



SHOW A RADIANT SMILE

WITH THE DENTACARE ORAL HYGIENE
SYSTEM FROM PANASONIC**Reliable cleaning power for your mouth**

Regular oral care is the most effective way to prevent tooth decay and periodontal disease. But even if you brush your teeth several times a day, your teeth can still be dangerously damaged or your gums become inflamed (gingivitis). This can even be caused by using the wrong cleaning technique or having a poor system of oral hygiene.

Be on the safe side with DentaCare sonic toothbrushes and oral irrigators from Panasonic. For a level of hygiene and cleanness you can see.

SONIC TECHNOLOGY FOR THE
ULTIMATE IN ORAL HYGIENE**Intensive cleaning with
Panasonic DentaCare sonic toothbrushes**

The process of cleaning your teeth with a manual toothbrush is quite individual. Brushing techniques, movement, power exerted and frequency of brush head movement can all vary, depending on the dexterity of the user and the length of time taken. As a result: even teeth that are cared for regularly can be brushed with the wrong technique, for too short a period of time, or with too much pressure. With DentaCare Sonodent, Panasonic has developed toothbrushes with ultra-modern sonic technology, to ensure optimum oral hygiene. At the same time, they protect and care for your gums. Quite automatically.

**SUPERFAST VIBRATIONS**

The DentaCare Sonodent motor moves the brush head up and down at up to 34,000 vibrations* per minute. This makes it around ten times faster than a normal electric toothbrush. The motor can maintain this rate, whatever the battery charge, or however intensively you brush.

**HORIZONTAL MOTION**

Because of their minimal vibration displacement amplitude (0.6-1mm), Panasonic sonic toothbrushes are very easily directed and positioned. The compact brush head can reach the spaces between your teeth that would otherwise be difficult to clean. Plaque and discoloration are easily removed.

**GENTLE BRUSH PRESSURE**

At 100-150g*, DentaCare Sonodent sonic toothbrushes exert only half as much pressure on your gums as many conventional electric toothbrushes. This avoids possible gum irritation and damage to the neck of the tooth.

*depending on the model

THE SONIC TOOTHBRUSH FOR BRILLIANT CLEANING

Type: Panasonic DentaCare Sonodent EW1031

Intensive cleaning power – at an unbelievably good price. At 26,000 vibrations per minute, the EW1031 DentaCare Sonodent with silver finish not only removes plaque completely, it adds glamour to your bathroom. Thanks to its new type of brush attachment with superfine microbristles, the sonic toothbrush also gets to places that were inaccessible to your old brush.

- Ultra-modern motor (26,000 vibrations/minute)
- Around 40 minutes of continuous operation
- Works at full power right up to recharging
- A non-contact charging station makes it easy to clean
- Includes a charging station with brush stand and two EW0911 brushes
- 2-year warranty



SONODENT EW1031 TECHNICAL DATA

Mains voltage: ~ 100-240 V; 50-60 Hz
Power: 2 W (fully charged)
Service life: approx. 40 min. continuous operation

DentaCare
Sonodent



PZN 1084789

THE SONIC TOOTHBRUSH WITH A POWERFUL LITHIUM ION BATTERY

Type: Panasonic DentaCare Sonodent EW-DL40

Get dazzlingly beautiful teeth and save the environment. The new lithium ion battery keeps the charging time short and the service life long. This saves energy. At the same time, you benefit from the intensive cleaning power of up to 31,000 vibrations/minute.

- Lithium ion battery: energy-saving short charging time and long service life
- Intensive and reliable cleaning at up to 31,000 vibrations/minute
- Two levels of intensity – normal and gentle – to individualize your oral care
- Excess pressure prevention function to protect your gums
- 2-minute brushing timer with 30-second intervals to maintain the recommended brushing time
- Includes carrying case, charging station with brush stand and an EW0925 precision-tip brush
- 2-year warranty

**NEW from
October
2009**



SONODENT EW-DL40 TECHNICAL DATA

Powerful lithium ion battery
Mains voltage: ~ 100-240 V; 50-60 Hz
Power: 2 W (fully charged)
Service life: approx. 60 min. continuous operation

DentaCare
Sonodent



PZN 00025052

THE SONIC TOOTHBRUSH WITH A MULTI-DIRECTIONAL BRUSH HEAD

Type: Panasonic DentaCare Sonodent EW1035

Cleaning your teeth has never been this versatile. The brushes of the EW1035 DentaCare Sonodent move sideways, horizontally and vertically at high speed, to ensure the most thorough cleaning, even in the most remote corners of your mouth. Different brush head versions are available, to assist with your oral hygiene. For an even cleaner and more radiant result.



- Ultra-modern linear motor (up to 34,000 vibrations/minute)
- Cleaning at its most thorough: diagonally, horizontally and vertically
- Circulation-enhancing massage mode
- Soft-start function
- Excess pressure warning function
- 2-minute timer
- Works at full power right up to recharging
- Includes six brush attachments, carrying case and storage box for brushes and charging station



SONODENT EW1035 TECHNICAL DATA

Mains voltage: ~ 100-240 V; 50-60 Hz
Power: 2 W (fully charged)
Service life: approx. 40 min. continuous operation

DentaCare
Sonodent




PZN 4780756

Optimum, multi-directional cleaning technique



With 31,000 horizontal and 34,000 vertical vibrations every minute, the EW1035's innovative linear motor removes tooth deposits gently and effectively. Precise, 1mm vibrations ensure that cleaning is effective yet gentle on your gums.



Brush and protect functions



Timer

The Panasonic timer* helps you keep to the brushing time dentists recommend. As a reminder, the toothbrush stops working briefly after 2 minutes and then starts up again. The EW-DL40 gives you additional help with 30-second intervals.



Soft-start

The soft-start function gives a more gentle start to brushing. It helps you protect your teeth and gums.

The vibration intensity automatically reduces if you press your toothbrush too hard against your teeth or gums. This is particularly kind to your gums.



All Panasonic toothbrushes and their chargers are fully washable, which makes them particularly easy to clean.

*depending on the model

DENTAL CARE WITH SMALL HEADS

Brush attachments to individualize oral hygiene

Not every tooth is upright and the spaces between your teeth also vary, depending on the nature of your mouth. Panasonic has seven different attachments to offer for your quite individual oral hygiene routine.

DentaCare brush attachments fit all Panasonic sonic toothbrushes.

The effective three-head toothbrush

Type: Panasonic Sonodent EW0900

Standard accessory for EW1012

- Thorough cleaning of all tooth surfaces at the same time
- Particularly recommended for back molars or for use by children and adults with restricted movement



PZN 1019504

The gentle, spring-loaded brush

Type: Panasonic Sonodent EW0920

Standard accessory for EW1012

- Intensive yet gentle cleaning from spring-loaded bristles
- The small brush head also lets you clean those difficult to reach interdental spaces



PZN 1019510

The brush with superfine microbristles

Type: Panasonic Sonodent EW0911

Standard accessory for EW1031

- The bristles have microfine ends for even better cleaning of interdental spaces



PZN 1822862

The precision-tip brush

Type: Panasonic Sonodent EW0925

Standard accessory for EW-DL40

- Longer bristles at the edge allow perfect adaptation to uneven tooth surfaces
- The spaces between your teeth are easily cleaned



PZN 0025049

The zig-zag bristle brush

Type: Panasonic Sonodent EW0923

Standard accessory for EW1035

- Microfine zig-zag bristles make it easy to clean those difficult to reach interdental spaces
- The spring-loaded bristles are gentle on your teeth and gums



PZN 4928526

The pin brush for in-depth cleaning

Type: Panasonic Sonodent EW0940

Standard accessory for EW1035

- The pin-shaped brush for thorough cleaning of deeper interdental spaces



PZN 4928503

The silicone brush for massaging gums

Type: Panasonic Sonodent EW0931

Standard accessory for EW1035

- The rubber-coated brush for a circulation-enhancing massage of your gums and the roof of your mouth



PZN 4928495



**NEW from
October
2009**

STATIONARY ORAL IRRIGATOR FOR USING WITH MOUTHWASH

Type: Panasonic DentaCare EW-DJ40

Is mouthwash part of your daily oral hygiene routine? Then the EW-DJ40 is exactly what you need! Just use your mouthwash* to clean your interdental spaces and enjoy the fresh, clean feeling it gives you. The leftover food between your teeth is easily and conveniently removed by a jet of water. You can adjust the jet intensity to suit your individual needs. The water tank is integrated in the handle and the charger is very small. Which means that the EW-DJ40 does not take up much space in your bathroom.

- Extremely easy to maneuver, it can be telescoped to a mere 14 cm
- 2 levels of water jet intensity
- 165 ml tank capacity (25% increase**)
- Pump frequency: 1,500 pulses per minute
- Reduced charging time – a mere 8 hours
- Includes a charger and an EW0955 nozzle

THOROUGH CLEANING OF YOUR INTERDENTAL SPACES

Stubborn leftover food in your interdental spaces. That unclean feeling in your mouth. You know what I mean? Then Panasonic oral irrigators are exactly what you need. DentaCare oral irrigators can clean the places your toothbrush cannot reach. Remove leftover food from between your teeth with a jet of water that is perceptibly thorough. This also makes them ideal for caring for your dentures or fixed braces. As well as thoroughly cleaning your teeth, they also gently massage and care for your gums.

Always use your oral irrigator immediately after brushing your teeth. You will be delighted by how fresh and clean your mouth will feel.

**NOW for using
with mouthwash!**

*diluted
**compared to EW1211

HANDY EW-DJ40 TECHNICAL DATA

Mains voltage: ~ 100-240 V; 50-60 Hz
Pump frequency: 1,500 pulses/min.
Service life:
approx. 15 min. continuous operation/
approx. 40 sec. per tank

DentaCare Handy



PZN 0025129

**NOW for
using with
mouthwash!**

STATIONARY ORAL IRRIGATOR

Type: Panasonic DentaCare EW1211

The DentaCare EW1211 oral irrigator is the compact addition to your oral hygiene routine. The water jet removes leftover food from areas that even the best toothbrush cannot reach and its gentle massage pressure enhances circulation in your gums. A second pump allows you to choose between a powerful jet of water and a massaging water-and-air mixture.

- Gentle cleaning with a dual air and water pump
- Water tank integrated in the handle
- Ergonomic design
- 3 levels of water jet intensity
- Pump frequency: 1,400 pulses per minute
- The 135 ml tank lasts around 35 seconds
- A non-contact charging station makes it easy to clean
- Includes charging station and two nozzles EW0950



HANDY EW1211 TECHNICAL DATA

Mains voltage:
~100-240 V; 50-60 Hz
Pump frequency: 1,400 pulses/min.
Service life:
approx. 15 min. continuous operation/
approx. 35 sec. per tank

DentaCare Handy



PZN 1025597

TRAVEL ORAL IRRIGATOR FOR USING WITH MOUTHWASH

Type: Panasonic DentaCare EW-DJ10

Now you can also use the new EW-DJ10 travel oral irrigator with mouthwash. You also benefit from a 20%* increase in water pressure – for even more thorough cleaning and that fresh feeling. The water tank of the EW-DJ10 is integrated in the handle. Which makes it very easy to maneuver. It is battery-operated, so you do not have to rely on a power socket. You never have to do without perfect oral and dental hygiene, even when travelling, or in a cramped bathroom.

- Extremely easy to maneuver, it can be telescoped to a mere 14 cm
- Battery-operated – no bulky power supply unit required
- 2 levels of water jet intensity
- Elliptical in shape, it is easy to hold and store
- Pump frequency: 1,400 pulses per minute
- The 165 ml tank lasts around 40 seconds
- Weighs only 210g
- Includes two batteries and an EW0955 nozzle

*compared to EW1270

TRAVEL EW-DJ10 TECHNICAL DATA

Power supply: 3 V DC voltage; 50-60 Hz;
two 1.5 V alkali batteries [AA/LR6]
Pump frequency: 1,400 pulses/min.
Service life:
approx. 20 min. continuous operation/
approx. 40 sec. per tank

DentaCare Travel



14 cm

PZN 0025112

SUPER-COMPACT TRAVEL ORAL IRRIGATOR

Type: Panasonic DentaCare EW1270

Cramped bathrooms or overfull suitcases are no longer a problem. The collapsible, battery-operated DentaCare Travel oral irrigator provides oral hygiene at home, on the move and on holiday. A tank integrated into the handle provides a lasting and invigorating water jet – as well as thorough cleaning.

- Extremely easy to maneuver, it can be telescoped to a mere 15 cm
- Battery-operated – no bulky power supply unit required
- Elliptical in shape, it is easy to hold and store
- Pump frequency: 1,600 pulses per minute
- The 165 ml tank lasts around 40 seconds
- Weighs only 210g
- Includes the battery and an EW0921 nozzle



TRAVEL EW1270 TECHNICAL DATA

Power supply: 3 V DC voltage; 50-60 Hz;
two 1.5 V alkali batteries [AA/LR6]
Pump frequency: 1,600 pulses/min.
Service life:
approx. 20 min. continuous operation/
approx. 40 sec. per tank



Nozzle attachments for thorough cleaning

For the spaces between your teeth to be cleaned effectively, the nozzles must work properly. The limescale in tap water may clog them up, over time. This will interfere with the efficiency of your oral irrigator. So we recommend that you replace your nozzles with new ones, about every six months.

Replacement nozzle for EW-DJ10 and EW-DJ40

Type: Panasonic EW0955

- 2 per pack
- Includes four rings in different colors. This makes it easier to tell the nozzles apart, if you share the use of your oral irrigator with several people



PZN 5467369

Replacement nozzle for EW1270

Type: Panasonic EW0921

- 2 per pack
- Colors: white and gray
- To make the nozzles more distinctive when the oral irrigator is being used by two people



PZN 1019533

Replacement nozzle for EW1211

Type: Panasonic EW0950

- 2 per pack
- Includes four different-colored rings, to make the nozzles more distinctive when the oral irrigator is being used by several people



PZN 0036452

THE DENTACARE ORAL HYGIENE SYSTEM FROM PANASONIC.

Our knowledge for your well-being.



Panasonic

Panasonic Deutschland

a division of Panasonic Marketing Europe GmbH

Winsbergring 15 · 22525 Hamburg, Germany

Telephone +49 (0) 40 - 85 49 0

Telefax +49 (0) 40 - 85 49 25 00

www.panasonic.de/mundpflege

Info hotline +49 (0) 180 - 51 51 800

(14 ct./min. using the network of Deutsche Telekom, status quo 09/2009)